


# Product data sheet



MedKoo Cat#: 462325 Name: NVS-PAK1-C CAS: 2250019-95-3 Chemical Formula: C <sub>23</sub> H <sub>25</sub> ClF <sub>3</sub> N <sub>5</sub> O Exact Mass: 479.17 Molecular Weight: 479.9322	
Product supplied as:	Powder
Purity (by HPLC):	≥ 98%
Shipping conditions	Ambient temperature
Storage conditions:	Powder: -20°C 3 years; 4°C 2 years. In solvent: -80°C 3 months; -20°C 2 weeks.

## 1. Product description:

NVS-PAK1-C is a potent, ATP-competitive and specific allosteric PAK1 inhibitor probe.

## 2. CoA, QC data, SDS, and handling instruction

SDS and handling instruction, CoA with copies of QC data (NMR, HPLC and MS analytical spectra) can be downloaded from the product web page under “QC And Documents” section. Note: copies of analytical spectra may not be available if the product is being supplied by MedKoo partners. Whether the product was made by MedKoo or provided by its partners, the quality is 100% guaranteed.

## 3. Solubility data

Solvent	Max Conc. mg/mL	Max Conc. mM
DMF	50.0	104.18
DMSO	50.0	104.18
DMSO:PBS (pH 7.2) (1:1)	0.5	1.04
Ethanol	20.0	41.67

## 4. Stock solution preparation table:

Concentration / Solvent Volume / Mass	1 mg	5 mg	10 mg
1 mM	2.08 mL	10.42 mL	20.84 mL
5 mM	0.42 mL	2.08 mL	4.17 mL
10 mM	0.21 mL	1.04 mL	2.08 mL
50 mM	0.04 mL	0.21 mL	0.42 mL

## 5. Molarity Calculator, Reconstitution Calculator, Dilution Calculator

Please refer the product web page under section of “Calculator”

## 6. Recommended literature which reported protocols for in vitro and in vivo study

In vitro study

TBD

In vivo study

TBD

## 7. Bioactivity

Biological target:

NVS-PAK1-C is a potent, ATP-competitive and specific allosteric PAK1 inhibitor probe with IC<sub>50</sub> values of 5 nM and 6 nM for dephosphorylated PAK1 and phosphorylated PAK1.

In vitro activity

TBD

In vivo activity

# Product data sheet



TBD

*Note: The information listed here was extracted from literature. MedKoo has not independently retested and confirmed the accuracy of these methods. Customer should use it just for a reference only.*